






# Inhibitor-Tolerant qPCR Mix

Optimized for diagnostic assays using clinical or environmental specimens with minimal sample pre-processing

## Suitable Application

-  Human Diagnostics
-  Vet Health
-  Environmental
-  Food Testing
-  Water Testing



A variety of different inhibitors inherent to clinical and environmental samples can negatively impact the sensitivity and accuracy of a molecular assay. Conventionally, to overcome this, specimens undergo processing prior to testing. However, extraction technologies are not 100% efficient which can impact the amount of target nucleic acid available for testing, and some inhibitors can co-elute following purification and induce false-negative assay results.

Inhibitor-resistant qPCR formulations are designed to reduce the amount of sample processing required for reliable qPCR analysis. Meridian Life Science has developed a novel inhibitor-resistant qPCR mix that can be used for amplification from crude lysates or inhibitor-rich samples such as urine, cerebral spinal fluid (CSF) or blood. It is suitable for developing qualitative assays that require fast turn-around times (TAT) or minimal sample processing to prevent sample loss.

## Inhibitor-Tolerant qPCR Mix


- Ready-to-use qPCR mix engineered with a proprietary buffer system designed for overcoming common PCR inhibitors
- Allows for DNA amplification direct from crude lysate or unprocessed sample such as urine, cerebral spinal fluid (CSF), milk and blood
- Ideal for developing qualitative assays where sample is limited
- Eliminates need for sample extraction and speeds up assay TAT

PRODUCT	CAT NO.	VOLUME	REACTIONS
Inhibitor-Tolerant qPCR Mix	MDX013	5 mL	500 Rxn
		100 mL	10,000 Rxn

## High qPCR efficiency across a range of inhibitors

Reaction efficiencies were determined from reactions containing a variety of known PCR inhibitors ranging from whole blood to biofluids (20% in reaction). The results demonstrate that the reaction efficiency of the Inhibitor-Tolerant qPCR Mix remained within 90-110% in the presence of a wide range of common PCR inhibitors.


**Inhibitor-Tolerant qPCR Mix Reaction Efficiencies**



## Ideal for developing quantitative assays direct from clinical or environmental samples


**Figure 1. Efficient amplification from various liquid sample types (spinal fluid, urine and milk).**

A 10-fold serial dilution of genomic DNA was spiked into cerebrospinal fluid, human urine or cow whole milk and the DNA was amplified, using Inhibitor-Tolerant qPCR Mix (green) and KAPA Probe Force Mix (red). The results illustrate that Inhibitor-Tolerant qPCR Mix is more sensitive than KAPA Probe Force Mix, as lower dilutions could be detected with better efficiencies, across all three sample types.



**Figure 2. Direct amplification from plant lysate.**

Genomic DNA was amplified from tomato lysate (1.6 mg tomato leaf tissue per reaction) in a triplex reaction. The results illustrate that Inhibitor-Tolerant qPCR Mix is able to amplify from samples rich in phenolic compounds, such as tomato crude leaf lysate, which are potent inhibitors of PCR.



### Ordering information:

**USA**  
E: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)  
Toll free: +1 800 327 6299

**Germany**  
E: [info.de@bioline.com](mailto:info.de@bioline.com)  
Tel: +49 (0)3371 60222 00

**美国迈迪安生命科学公司**  
电子邮件: [vivian.li@meridianlifescience.com](mailto:vivian.li@meridianlifescience.com)  
电话: +86-159-1103-0750

**UK**  
E: [info.uk@bioline.com](mailto:info.uk@bioline.com)  
Tel: +44 (0)20 8830 5300

**France**  
E: [info.fr@bioline.com](mailto:info.fr@bioline.com)  
Tel: +33 (0)1 42 56 04 40

**Australia**  
E: [info.au@bioline.com](mailto:info.au@bioline.com)  
Tel: +61 (0)2 9209 4180