

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: 9552 **Lot #:** 6G19120

Description: SARS-CoV-2 Spike S1 (RBD)
SARS-CoV-2 Spike S1 RBD (Receptor-Binding Domain) Protein, Recombinant
His-Tag at C-Terminus

Source: HEK293

Format: Affinity Purified, Liquid

Purity: ≥ 90% (SDS-PAGE)

Concentration: 2.00 mg/mL (BCA)

Buffer: 1X Phosphate Buffered Saline, pH 7.2

Preservative: None

Applications: Suitable for use in ELISA, Lateral Flow and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store unopened at -80°C. Working aliquots can be prepared and stored at or below -20°C for three months. Avoid repeated freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Quality Signature:



09 JUL 2020

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY