

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: 9559 **Lot #:** 3F15420

Description: SARS-CoV-2 Spike S2 (His-Tag)
SARS-CoV-2 Spike Protein, S2 Subunit (His-Tag), Recombinant

Source: Insect Cells

Format: Affinity Purified

Purity: ≥ 85% (SDS-PAGE)

Concentration: 1.9 mg/mL (BCA)

Buffer: Phosphate Buffered Saline, pH 7.4

Preservative: None

Applications: Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at -20°C. Avoid repeated freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Quality Signature:



02 JUN 2020

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY