



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01409H **Lot #:** 36L33815

Description: 1,25-Dihydroxy Vitamin D3

Source: Synthetic

Format: Purified, Lyophilized
Reconstitute in 99% Ethanol.

Purification: > 99% pure (HPLC).

Concentration: Not determined

Buffer: Lyophilized from an Ethanol Solution.

Preservative: None

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at -20°C protected from heat and light.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Sincerely, J. Orr

Signature

11 July 2017

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY