



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** A01429H **Lot #:** 3G19114

**Description:** Troponin I-C Complex, Recomb.  
Troponin I Cardiac (Fragment) and Troponin C Cardiac (Full Length) Complex, Recombinant  
Contains 6 Histidine Tag

**Source:** *Pichia pastoris*

**Format:** Purified, Liquid

**Purification:** > 85% pure (determined by SDS-PAGE).

**Concentration:** 1.3 mg/mL (Determined by BCA.)

**Buffer:** PBS, pH 7.4

**Preservative:** 0.1% Sodium Azide

**Applications:** Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Aliquot and store at -70°C. Avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Signature

13 May 2016

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

002