



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: A01720H **Lot #:** 4K30517

Description: Human Myeloma IgE lambda
Human Myeloma Immunoglobulin E (IgE) lambda

Source: Normal Human Plasma spiked with cell culture Human IgE, lambda.

Format: Neat, Liquid

Purification: Human Myeloma IgE lambda in neat Human plasma (EDTA).

Concentration: 43.8 µg IgE/mL by Siemens Centaur CP which is standardized against WHO 75/502.
Activity: 18,269 IU/mL (1 IU = 2.4 ng IgE) (Centaur CP)

Buffer: Not Applicable

Preservative: None

Applications: Specific applications have not been tested using this product, but use in such applications should not necessarily be excluded.

Storage: Upon receipt, aliquot and store below -10°C. Avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Brenda Dunn

01Nov2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY