



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** A50796H **Lot #:** 2D11916

**Description:** Human Myeloma IgE, lambda  
Human Myeloma Immunoglobulin E (IgE), lambda

**Source:** Human Myeloma Plasma

**Format:** Purified, Liquid

**Purification:** > 95% pure (SDS-PAGE).  
Single arc by IEP against antisera to whole human serum, human IgE and human lambda light chain.

**Concentration:** 1.18 mg/mL

**Buffer:** 15 mM Potassium Phosphate Buffer, pH 7.4 containing 150 mM NaCl.

**Preservative:** 0.05% Sodium Azide

**Applications:** Suitable for use in Immunoelectrophoresis. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Upon receipt, aliquot and store at  $\leq -70^{\circ}\text{C}$ . Avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*CORRY THOMAS*

22 Jun 2017

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

002