

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: B65510R **Lot #:** 1C06821

Description: Rabbit anti Group A Strep.
Rabbit Antibody to *Streptococcus* Group A

Specificity: Type specific carbohydrate of group A *Streptococcus*. Does not react with other Strep groups.

Host Animal: Rabbit

Immunogen: *Streptococci*, Group A

Format: Affinity Purified, Liquid

Purification: Affinity Chromatography

Concentration: 1.0 mg/mL (OD280nm, E^{0.1%} = 1.4)

Buffer: 0.01 M Phosphate Buffered Saline, pH 7.2
No stabilizing proteins have been added.

Preservative: 0.1% Sodium Azide

Applications: Suitable for use in ELISA, IFA and Lateral Flow immunoassays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pair for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
C01756R	B65510R
B65510R	C01756R

Storage: Short-term (up to 6 months) store at 2–8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Quality Signature:



10 MAR 2021

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY