



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: C01247M **Lot #:** 6F17018

Description: MAb to HBsAg (ad/ay)
Monoclonal Antibody to Hepatitis B Surface Antigen (HBsAg)

Specificity: Recognizes Hepatitis B surface antigen, subtypes ad & ay. Does not cross react with normal human serum.

Host Animal: Mouse **Isotype:** IgG_{2a}

Source: Ascites

Immunogen: Purified HBsAg

Format: Purified, Liquid

Purification: > 90% pure. Protein A Chromatography

Concentration: 3.3 mg/mL (OD280nm, E^{0.1%} = 1.3)

Buffer: 0.01 M Phosphate Buffered Saline, pH 7.2
This product contains no stabilizing protein.

Preservative: 0.1% Sodium Azide

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pairs for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
C01247M	C01249M
C01247M	B65811P

Storage: Store at 2-8°C.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Brenda Dunn

21 June 2018

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY