



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road

Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: C01873M **Lot #:** 3D11617

Description: MAb to Norovirus GI Capsid
Monoclonal Antibody to Norovirus GI Capsid Protein

Specificity: Recognizes the GI capsid protein of Norovirus (native and recombinant). Does not cross react with Norovirus GII Capsid Protein.

Host Animal: Mouse **Isotype:** IgG_{2a}

Source: Ascites

Immunogen: Recombinant Norovirus GI Capsid Protein

Format: Purified, Liquid

Purification: > 90% pure. Protein A Chromatography

Concentration: 2.5 mg/mL (OD280nm, E^{0.1%} = 1.3)

Buffer: 0.01 M Phosphate Buffered Saline, pH 7.2
Product contains no stabilizing proteins.

Preservative: 0.1% Sodium Azide

Applications: Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pairs for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
C01873M	C01875M

Storage: Store at 2-8°C.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

CORRY THOMAS

26 April 2017

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY