



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** C01876M **Lot #:** 2G20916

**Description:** MAb to Norovirus GII Capsid  
Monoclonal Antibody to Norovirus GII Capsid Protein

**Specificity:** Recognizes the GII capsid protein of Norovirus (native and recombinant). Does not cross react with Norovirus GI Capsid Protein

**Host Animal:** Mouse **Isotype:** IgG<sub>1</sub>

**Source:** Ascites

**Immunogen:** Recombinant Norovirus GII Capsid Protein

**Format:** Purified, Liquid

**Purification:** > 90% pure. Protein A Chromatography

**Concentration:** 3.1 mg/mL (OD280nm, E<sup>0.1%</sup> = 1.3)

**Buffer:** 0.01 M Phosphate Buffered Saline, pH 7.2  
Product contains no stabilizing proteins.

**Preservative:** 0.1% Sodium Azide

**Applications:** Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pairs for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
C01876M	C01877M

**Storage:** Store at 2-8°C.

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Brenda Dinn*

11Dec2017

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**