



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: C01903M **Lot #:** 3C06817

Description: MAb to Norovirus
Monoclonal Antibody to Norovirus

Specificity: Recognizes human norovirus, feline calicivirus and rabbit hemorrhagic disease virus.

Host Animal: Mouse **Isotype:** IgG_{2b}

Source: Ascites

Immunogen: Rabbit hemorrhagic disease virus and feline calicivirus.

Format: Purified, Liquid

Purification: > 90% pure. Protein A Chromatography

Concentration: 3.1 mg/mL (OD280nm, E^{0.1%} = 1.4)

Buffer: Phosphate Buffered Saline, pH 7.4

Preservative: 0.09% Sodium Azide

Applications: Suitable for use in ELISA, Lateral Flow, Hemagglutination Inhibition Assay and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pairs for immunoassays:

<u>Capture</u>	<u>Detection</u>	
C01902M	C01903M	Lateral Flow and ELISA
C01903M	C01904M	Lateral Flow and ELISA

Storage: Store at 2-8°C.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Signature

09 March 2017

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY