



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: C86220M **Lot #:** 2H23119

Description: MAb to Staph. aureus Enterotoxin B
Monoclonal Antibody to Staphylococcus aureus Enterotoxin B

Specificity: Staphylococcus aureus Enterotoxin B. Does not cross-react with Staphylococcus aureus enterotoxins A, C, D and E.

Host Animal: Mouse **Isotype:** IgG₁

Source: Ascites

Immunogen: Recombinant Staphylococcus aureus enterotoxin B expressed in E. coli.

Format: Purified, Liquid

Purification: ≥ 95% pure (SDS-PAGE). Protein G Chromatography

Concentration: 4.8 mg/mL (OD280nm, E^{0.1%} = 1.4)

Buffer: Phosphate Buffered Saline, pH 7.4

Preservative: 0.09% Sodium Azide

Applications: Detection of staphylococcal enterotoxins in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody pair for immunoassay of St. aureus enterotoxin B.

Capture

C86220M

Detection

C86400M

Storage: Store at 2-8°C.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Brenda Dinn

20AUG2019

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY