

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: E01310M **Lot #:** 4A01121

Description: MAb to Thyroid Peroxidase
Monoclonal Antibody to Human Thyroid Peroxidase (TPO)

Specificity: Human Thyroid peroxidase, native and recombinant.

Host Animal: Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from BALB/c mice. **Isotype:** IgG₁

Source: Ascites

Immunogen: Human Thyroid Peroxidase

Format: Purified, Liquid

Purification: > 90% pure (SDS-PAGE). Protein A Chromatography

Concentration: 6.3 mg/mL (OD280nm, E^{0.1%} = 1.4)

Buffer: Phosphate Buffered Saline, pH 7.4

Preservative: 0.09% Sodium Azide

Applications: Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

<u>Capture</u>	<u>Detection</u>
E01310M	E01309M

Storage: Store at 2-8°C.

Safety Notes (s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Quality Signature:



11 JAN 2021

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY