



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: G01306M **Lot #:** 8I26114

Description: MAb to EDDP
Monoclonal Antibody to 2-Ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine (EDDP)

Specificity: EDDP

Host Animal: Mouse **Isotype:** IgG₁

Immunogen: EDDP conjugated to Bovine Serum Albumin (BSA).

Source: Ascites

Format: Purified, Liquid

Purification: ≥ 95% pure (SDS-PAGE). Protein G Chromatography

Concentration: 5.23 mg/mL (OD280nm, E^{1%} = 1.4)

Buffer: Phosphate Buffered Saline, pH 7.4

Preservative: 0.02% Sodium Azide

Applications: Suitable for use in ELISA and Lateral Flow assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody-antigen pair for ELISA and Lateral Flow assays:

<u>Antigen-Coating</u>	<u>Antibody-Label</u>
Y01305B	G01306M
Y01306B	G01306M

Storage: Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Signature

09 February 2017

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY