



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: H8P01-767 **Lot #:** 4A02916

Description: HSA (Low B12/Low Folate)
Human Serum Albumin (HSA), Low B12 and Low Folate

Source: Phoresis Plasma

Format: Purified, Lyophilized

Purification: 96.4% pure (SPE gel electrophoresis). Cohn-Oncley Cold Alcohol Procedure

Concentration: Total Protein: 4.6 g/dL (Biuret)

pH: 7.16 (5% solution at room temperature).

Buffer: Not Applicable

Preservative: None

Applications: Specific methodologies have not been tested using this product.

QC Information: B-12: 40 pg/mL (Centaur)
Folate: BDL ng/mL (Centaur)
Bioburden (5%): < 100 cfu/mL by Pour Plate at 30-35°C after 48 hours.

Storage: Short term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Sincerely,
S. Orr*

04 January 2019

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

002