



Meridian

Life Science, Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
Telephone: 901-382-8716
Fax: 901-333-8223
Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: K24980R **Lot #:** 15D11303

Description: Rabbit anti Resistin (33-90)
Rabbit Antibody to Human Resistin (amino acids 33-90)

Specificity: Human Resistin

Host Animal: Rabbit

Immunogen: Synthetic Human Resistin peptide (a.a. 33-90) bTG conjugated
(CQSVTSRGDLATCPRGFAVTGCTCGSACGSWDVRAETTCHCQCAGMDWTGARCCRVQP)

Format: Neat, Lyophilized.
Reconstitute in 20ul double distilled water.

Purification: Not applicable. **Exp. Date:** N/A

Concentration: Not applicable.

Buffer: Not applicable.

Preservative: None

Applications: Suitable for use in ELISA.

Storage: Store lyophilized product at 2-8°C. After reconstitution, store at -20°C. Avoid multiple freeze/thaw cycles.

Robert Ott

Signature

23 Jan 2015

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY