



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: M01326M **Lot #:** 2I24819

Description: MAb to Cyfra-21-1
Monoclonal Antibody to Cyfra-21-1

Specificity: Reacts with Cyfra-21-1

Source: Mouse **Isotype:** IgG₁

Immunogen: Purified Human Cyfra-21-1 derived from tissue culture.

Format: Purified, Liquid

Purification: Ion Exchange Chromatography

Concentration: 2.7 mg/mL

Buffer: Phosphate Buffered Saline, pH 7.5

Preservative: 0.05% Sodium Azide

Applications: Suitable for use in ELISA applications. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody pair for sandwich assay:

<u>Capture</u>	<u>Detection</u>
M01326M	M01332M

Storage: Store at 2-8°C.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Brenda Dinn

05SEPT2019

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY