



Meridian

Life Science, Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716
Fax: 901-333-8223

Email: info@meridianlifescience.com
www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: MAG40-368 **Lot #:** H0950

Description: MAb to IgE
Monoclonal Antibody to Human Immunoglobulin E (IgE)

Specificity: Cross Reactivity:
Human IgG < 0.0001%
Human IgM < 0.0001%

Host Animal: Mouse **Isotype:** IgG₁, kappa

Source: Ascites

Format: Purified, Liquid **Exp. Date:** 08 NOV 2021

Purification: > 90% pure (SDS-PAGE). Protein A Chromatography

Concentration: 5.2 mg/mL (OD280nm, E^{0.1%} = 1.4)

Buffer: 10 mM Phosphate, pH 7.4 containing 150 mM Sodium Chloride.

Preservative: 0.1% Sodium Azide


Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody pair for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
MAG40-231	MAG40-368

Storage: Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.


Signature

08 November 2019
Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY