

Lyo-Compatible MMLV-RT

Suitable for Research and further Manufacturing Use

Catalog No:	MDX042
Lot No:	EN055-B092750
Storage Conditions:	-20°C
Component Lot No.	LCR-021302B
Expiry date:	March 2023

Quality Control Parameters

High-concentration MMLV-RT suitable for incorporation into lyophilized RT-PCR assays

Analysis	Specification	Result
Functional	<p>Activity is measured as reverse transcriptase units by quantitative PCR analysis against a reference enzyme.</p> <p><u>Pass Criteria:</u> Activity must be greater than 165 U/μL</p>	591 U/μL
DNA contamination	<p>Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in concordance with control sample.</p> <p><u>Pass Criteria:</u> Amplification traces must overlay with the negative control.</p>	Passed
DNase contamination	<p>DNase contamination is measured as DNA substrate degradation against a DNase I dilution series by agarose gel electrophoresis.</p> <p>Limit of detection: 6.25×10^{-4} KU DNase I.</p> <p><u>Pass Criteria:</u> No detectable degradation.</p>	Passed

United Kingdom

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

Australia

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763

RNase contamination	Quantitative PCR analysis with high and low RNase standards. Limit of detection: 9.7×10^{-3} ng/ μ L RNase <u>Pass Criteria:</u> Test sample must show less RNase activity than the limit of detection.	Passed
---------------------	---	--------

QA / QC Representative:



Andrew Galeeba-M

Date: 17th February 2021

United Kingdom

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

Australia

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763

Enzyme Dilution Buffer

For research or further manufacturing use only

Catalog No:	MDX042
Lot No:	EN055-B092750
Storage Conditions:	-20°C
Lot number:	TDB-021202B
Expiry date:	March 2023

Quality Control Parameters

Enzyme Dilution Buffer is a glycerol-free 10x buffer used for the dilution of Taq antibody or reverse transcriptase.

Analysis	Specification	Result
Functional	<p>A 3Kb fragment was amplified with a dilution series of Enzyme Dilution Buffer, using standard conditions and 30 cycles.</p> <p>A 3Kb fragment was amplified with a dilution series of <i>Taq polymerase</i>, using standard conditions and 30 cycles.</p> <p>Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).</p>	Passed
Glycerol determination	Glycerol content is < 0.2% determined by spectrophotometric analysis and comparison to a standard curve.	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 6.25×10^{-4} KU/ μ L.	Passed
RNase contamination	Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection: 9.7×10^{-3} ng/ μ L RNase.	Passed

QA / QC Representative:



Andrew Galeeba-M

 Date: 17th February 2021

United Kingdom

 Tel: +44 (0)20 8830 5300
 Fax: +44 (0)20 8452 2822

USA

 Tel: +1 901.382.8716
 Fax: +1 901.382.0027

Germany

 Tel: +49 (0)3371 60222 00
 Fax: +49 (0)3371 60222 01

Australia

 Tel: +61 (0)2 9209 4180
 Fax: +61 (0)2 9209 4763