



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road

Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** N62934H **Lot #:** 3I24615

**Description:** Apo AI/B Free Serum  
Apolipoprotein AI/B Free Serum

**Source:** Pooled Defibrinated Human Serum

**Format:** Purified, Liquid

**Purification:** Extraction of Apo AI and Apo B from normal human serum.  
The resulting serum is non-measurable for Apo AI and Apo B by TIA (Diasorin test kit).  
Product is 0.2 µm filtered.

**Concentration:** Total Protein: 4.2 g/dL (Biuret: Cobas-Mira)

**Buffer:** Not Applicable

**Preservative:** 0.1% Sodium Azide

**Applications:** Specific methodologies have not been tested using this product.

**Storage:** Short term store at 2–8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Warnings:** All human source materials have tested negative for HIV and HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.

This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.

Signature

03 SEP 2015

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**