



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road

Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** N86710M **Lot #:** 3I26115

**Description:** MAb to Fibrinopeptide A  
Monoclonal Antibody to Fibrinopeptide A

**Specificity:** Fibrinopeptide A region in the alpha-chain of fibrinogen and free fibrinopeptide A.

**Host Animal:** Mouse **Isotype:** IgG<sub>2a</sub>

**Source:** Ascites

**Immunogen:** Human Fibrinopeptide A, conjugated with a carrier protein.

**Format:** Purified, Liquid

**Purification:** > 90% pure (SDS-PAGE). Protein A Chromatography

**Concentration:** 7.3 mg/mL (OD280nm, E 0.1% = 1.4)

**Buffer:** Phosphate Buffered Saline, pH 7.4

**Preservative:** 0.09% Sodium Azide

**Applications:** Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pair for detection of fibrinogen in sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
N86710M	N86413M

**Storage:** Store at 2-8°C.

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Brenda Dinn*

16Oct2017

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**