



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: R01712 **Lot #:** 5A02220

Description: AMH Recombinant
Anti Mullerian Hormone Antigen, Recombinant
Contains a His-tag. Molecular weight ~32 kDa (partial AMH sequence).

Source: *E. coli*

Format: Purified, Liquid

Purification: > 90% pure (SDS-PAGE).

Concentration: 2.62 mg/mL

Buffer: 50 mM Tris, pH 8.0 with 0.15 M Sodium Chloride and 10% Glycerol.

Preservative: 0.02% Proclin 300

Applications: Suitable for use in sandwich and competitive ELISA, CLIA, Lateral Flow and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended Antibodies and Antigen for sandwich immunoassay:

<u>Capture MAb</u>	<u>Detection MAb</u>	<u>Recombinant Antigen</u>
E01347M	E01348M	R01712

Storage: Upon receipt, store at 2-8°C up to 1 week and at -20°C long term.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Blanca Gelle

22 JAN 2020

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY