



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: R01714 **Lot #:** 8H22618

Description: RSV, Long Strain
Respiratory Syncytial Virus (RSV) Antigen, Long Strain

Source: Vero Cell Culture
Strain: Long Strain (ATCC VR-26)

Format: Purified, Liquid

Purification: Chemical Extraction and Sonication

Concentration: 0.94 mg/mL (Bradford method using BSA as standard).

Buffer: 100 mM Glycine Buffer, pH 9.5

Preservative: None

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Store at $\leq -65^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles.

Inactivation: Beta-propiolactone treatment. Confirmation by attempted growth under original culture conditions.

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Dwanda Duman

17 Sept 2018

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY