



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** R01720 **Lot #:** 27K33217

**Description:** *C. coli* OMP Recombinant  
*Campylobacter coli* Outer Membrane Protein (OMP) Recombinant  
Product contains fusion partner His tag.  
Molecular Weight: 18 kDa

**Source:** *E. coli*

**Format:** Purified, Liquid

**Purification:** Multi-step procedure including affinity chromatography.

**Concentration:** 2.7 mg/mL (OD280nm)

**Buffer:** Tris Buffer

**Preservative:** 0.09% Sodium Azide

**Applications:** Suitable for use as control in Lateral Flow Assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

This antigen works with the following antibody pair:

Capture MAb      Detection MAb  
C01947M          C01946M

**Storage:** Store at 2-8°C for < 1 week. Long term storage at -20°C. Aliquot to avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Brenda Dinn*

29Nov2017

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**