



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: R18292 **Lot #:** 4B05516

Description: Rubella Virus E2, Recomb.
Rubella Virus E2 Protein, Recombinant

Contains the Rubella virus E2 region. Does not contain a fusion partner. Molecular weight is 33 kDa. Immunoreactive with Rubella virus positive sera.

Source: *E. coli*

Format: Purified, Liquid

Purification: > 95% pure (10% PAGE, Coomassie Staining). Ion Exchange Chromatography

Concentration: 1 mg/mL

Buffer: 20 mM MES, pH 6.5 containing 8 M Urea, 100 mM, Sodium Chloride.

Preservative: None

Applications: Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage: Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Sincerely,
S. Orr*

Signature

26 April 2017

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

002