



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** Y01268B **Lot #:** 1A03018

**Description:** THC BSA Conjugate  
Tetrahydrocannabinol (THC) Bovine Serum Albumin Conjugate  
9-Tetrahydrocannabinol (THC) Bovine Serum Albumin Conjugate

**Source:** Not Applicable

**Format:** BSA, Liquid

**Concentration:** 13.24 mg/mL

**Buffer:** 10 mM Phosphate Buffered Saline, pH 7.4

**Preservative:** 0.1% Sodium Azide

**Applications:** Suitable for use in ELISA and Lateral Flow assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody-antigen pair for ELISA and Lateral Flow assays:

<u>Antibody</u>	<u>Antigen</u>
G01286M	Y01268B

**Storage:** Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

**Inactivation:** Not Applicable

**Safety Note(s):** Refer to the appropriate Safety Data Sheet (SDS) for additional information.

*Brenda Dunn*

30Jan2018

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**