



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: Y01306B **Lot #:** 4I27314

Description: EDDP BSA Conjugate
2-Ethylidene-1, 5-Dimethyl-3, 3-Diphenylpyrrolidine (EDDP) Bovine Serum Albumin (BSA) Conjugate

Source: Not Applicable

Format: BSA, Liquid

Concentration: 2.22 mg/mL

Buffer: PBS, pH 7.4

Preservative: 0.02% Sodium Azide

Applications: Suitable for use in ELISA and Lateral Flow assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antigen-antibody pair for Competitive Inhibition Immunoassays:

<u>Antigen-Coating</u>	<u>Antibody-Label</u>
Y01306B	G01306M

Storage: Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

Inactivation: Not Applicable

Safety Note(s): Refer to the appropriate Safety Data Sheet (SDS) for additional information.

Signature

27 May 2016

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

001