



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Catalog #:** D10410S **Lot #:** 1H21515

**Description:** Sheep anti Thyroxine (T4)  
Sheep Antibody to Thyroxine (T4)

**Specificity:** Thyroxine (T4)

**Host Animal:** Sheep

**Immunogen:** T4/Bovine Thyroglobulin

**Format:** Neat, Liquid

**Purification:** Not Applicable

**Concentration:** Not Determined

**Buffer:** Not Applicable

**Preservative:** 0.1% Sodium Azide

**Applications:** Specific methodologies have not been tested using this product.

**QC Information:** Not Applicable

**Storage:** Short-term store at 2-8°C. Long term store at -20°C.

**Warnings:** This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.

Signature

06 Aug 2015

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**