



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: info@meridianlifescience.com

www.MeridianLifeScience.com

CERTIFICATE OF ANALYSIS

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: F30369S **Lot #:** 189-22401

Description: Sheep anti Vancomycin
Sheep Antibody to Vancomycin

Specificity: Vancomycin. Minimal cross-reactivity with related compounds.

Host Animal: Sheep

Immunogen: Vancomycin-KLH Conjugate

Format: Neat, Liquid

Purification: Not Applicable

Concentration: Not Determined
Titer: 1:7,000

Buffer: Not Applicable

Preservative: None

Applications: Specific methodologies have not been tested using this product.

Storage: Upon receipt, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

Signature

19 September 2013

Date

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY