



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
Telephone: 901-382-8716  
Fax: 901-333-8223  
Email: info@meridianlifescience.com  
www.MeridianLifeScience.com

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** K95002S **Lot #:** 13I2589

**Description:** Sheep anti Thrombin  
Sheep Antibody to Human Thrombin

**Specificity:** Reacts with Purified Thrombin and Prothrombin. Cross-reacts with Human and mouse.

**Host Animal:** Sheep

**Immunogen:** Purified Human Thrombin (derived from purified prothrombin).

**Format:** Purified, Liquid

**Purification:** Caprylic Acid Precipitation

**Concentration:** 18 mg/mL (OD280nm, E<sup>1%</sup> = 14)

**Buffer:** 10 mM HEPES, 0.05 M Sodium Chloride, pH 7.4 containing 50% Glycerol.

**Preservative:** None

**Applications:** Suitable for use in ELISA. May be used for affinity purification. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Upon receipt, store at -20°C. Avoid multiple/freeze thaw cycles.

*Robert Ott*

Signature

30 September 2013

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY.**