



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Catalog #:** R18401 **Lot #:** 3I27304

**Description:** VZV Ag, ORF9 (Dumas) Recomb.  
Varicella-Zoster Virus (VZV) ORF9, (Dumas Strain), Recombinant  
  
Contains the VZV ORF9 immunodominant regions. Contains GST fusion partner. Immunoreactive with VZV positive sera.

**Source:** *E. coli*  
Strain: Dumas

**Format:** Purified, Liquid

**Purification:** > 95% pure (10% PAGE Coomassie Staining). GS-4B Sepharose Affinity Purification.

**Concentration:** 1 mg/mL

**Buffer:** 25 mM Tris-HCl, 1 mM EDTA containing 50% Glycerol.

**Preservative:** None

**Applications:** Suitable for use in ELISA and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Upon receipt, store at -20°C.

**Inactivation:** Not Applicable

Signature

20 September 2013

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**