



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134

USA

Telephone: 901-382-8716

Fax: 901-333-8223

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)

[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

### CERTIFICATE OF ANALYSIS

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** W01238S **Lot #:** 16D10708.A

**Description:** Sheep anti Urokinase  
Sheep Antibody to Human Urokinase (uPA)

**Specificity:** Recognizes Human uPA both alone and in complex under reducing and native conditions.

**Host Animal:** Sheep

**Immunogen:** Native Human uPA

**Format:** Purified, Liquid

**Purification:** Protein G Chromatography

**Concentration:** 1 mg/mL (OD280nm,  $E^{0.1\%} = 1.36$ )

**Buffer:** 0.05 M Sodium Phosphate, pH 6.6 containing 0.1 M Sodium Chloride and 1 mM EDTA.

**Preservative:** None

**Applications:** Suitable for use in ELISA (1:100,000) and Western Blot (1-5 µg/mL). Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage:** Upon receipt, store at  $-70^{\circ}\text{C}$ . Avoid multiple freeze/thaw cycles.

Signature

18 September 2013

Date

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**