



Redefining the possibilities of life from discovery to diagnosis

40+ Years History

650+ Employees

\$200M+ in sales

12 Worldwide Locations

70+ Global Markets

Nasdaq: VIVO



Global locations

POINT-OF-CARE DIAGNOSTICS



Blood Lead Testing

Results in 3 min
Easy to use, portable finger prick test

Anodic stripping voltammetry technology

3 Steps to a quantitative result

- 1. Collect:** Fingerstick sample (50µL)
- 2. Mix:** Dispense into reagent
- 3. Test:** Insert sensor, apply sample



Rapid Lateral Flow

Results in under 15 min
Single use

Testing Menu

- Tropical Disease:** Malaria, Dengue, CHIKV, Zika, LEISH, LEPTO
- Cardiac:** cTNI, ck-MB, MYO, D-Dimer, FABP
- Inflammation:** CRP, PCT
- Cancer:** CPT, FOB
- Infectious Disease:** HAV, HEV, HIV, Syphilis, TB, *S. typhi*, *H. pylori*, *C. diff*, Noro, Rota/Adeno
- Woman's Health:** TORCH, Rubella, Toxo, CMV, HSV, CT/GC, GBS
- Respiratory Disease:** Flu, *M. pne.*, Adeno, Legionella, GAS, RSV



Molecular

Results in under 60 min
Flexible, 1-10 tests/run

Testing Menu

- C. difficile*
- M. pneumoniae*
- Group A Strep
- Group B Strep
- Pertussis
- HSV 1 & 2
- Malaria
- Chlamydia (CT)
- Gonorrhea (NG)



LIFE SCIENCE RAW MATERIALS



For human, animal and plant diagnostic assays

MOLECULAR

- PCR enzymes for end-point PCR, qPCR, RT-qPCR & isothermal amplification
- Ultra-pure nucleotides for commercial scale production
- Manufacturing in Germany

Key Enzymes

- **MyTaq™ HS DNA Polymerase:** Antibody mediated hot-start DNA polymerase
- **Izopol Thermostable DNA polymerase:** For isothermal amplification
- **SensiFAST™ Probe:** qPCR for cDNA/DNA template
- **SensiFAST™ Probe One-Step Kit:** RT-qPCR for RNA template
- **IMMOLASE™:** Chemical mediated hot-start DNA polymerase
- **Jet Seq™:** Next Generation Sequencing

Glycerol-free, Lyo-ready formulations available

IMMUNODIAGNOSTICS

- Infectious disease expertise
- +3400 antigens & antibodies
- Large lot sizes & lot reservations

Key Diseases

- **Infectious Disease**

- ToRCH
- Childhood
- STDs
- Gastrointestinal
- Food & Water
- Respiratory
- Tropical

- **Cancer**

- **Cardiac**

- **Hormones**

- **Allergens**

- **Vet Health**

- **Autoimmune**

- **Drugs of Abuse**

- **Blockers & Immunoglobulins**

Key Products:

TRU Block™ | Zika
Vitamin D | Rubella
Nipah | PIVKA-II
Mouse IgG
Dengue



To place an order, please contact:

www.MeridianBioscience.com/LifeScience

Meridian Life Science Inc.

5171 Wilfong Road
Memphis, TN 38134

901-382-8716 / 800-327-6299
info@meridianlifescience.com