

# dNTPs (Deoxynucleotide triphosphates)

Ideal for PCR, qPCR, RT-PCR, cDNA synthesis, NGS, isothermal amplification, hi-fidelity and long-range PCR.



## MOST ECONOMIC OPTION FOR COMMERCIAL SCALE

- Manufactured in hundreds of litres per year
- Value offer for large quantities



## FOR DIAGNOSTIC & APPLIED MARKETS

- Quality specifications for highly sensitive assays
- Sodium salt form, >99% (HPLC)
- Free from DNase, RNase and PCR inhibitors



## RELIABLE SUPPLY & READY TO SHIP

- Large quantities in stock
- Pre-shipment samples available
- Stable for 3 years at -20°C, >50 freeze thaw cycles

Cat #	Product Name	Cat #	Product Name
<b>MDX067</b>	dNTP Mix, 10mM sodium salt	<b>MDX064</b>	dCTP, 100mM sodium salt
<b>MDX083</b>	dNTP Mix, 40mM sodium salt	<b>MDX065</b>	dGTP, 100mM sodium salt
<b>MDX084</b>	dNTP Mix, 100mM sodium salt	<b>MDX066</b>	dTTP, 100mM sodium salt
<b>MDX063</b>	dATP, 100mM sodium salt	<b>MDX096</b>	dUTP, 100mM sodium salt

## Ordering information:

**USA**  
5171 Wilfong Road  
Memphis, Tennessee 38134

Fax: +1 901-333-8223  
Toll Free: +1 800 327 6299

Email: [info@meridianlifescience.com](mailto:info@meridianlifescience.com)  
**Orders: [orders@meridianlifescience.com](mailto:orders@meridianlifescience.com)**  
[www.MeridianLifeScience.com](http://www.MeridianLifeScience.com)

Connect with us:

ISO 13485 Certified

03/20